

Connecting your hearing aids to the Signia app for remote control access:

1. Begin by downloading the Signia app and tapping on it. You can access this from your app store (iphone), or play store (android).
2. Make sure your hearing aids are connected to Bluetooth and that Bluetooth is on. When you open the app, select the option for “Bluetooth.”
3. To enable the pairing mode, switch your hearing aids off and then back on. If you use disposable batteries, open the battery door and then close it. If you have rechargeable hearing aids, place in the charger for a few seconds and then pull them out. The instructions on the screen will prompt you through each step.
4. Once the app finds your aids, tap on the screen to select your aids.
5. After a few seconds the aids will be connected to the app, and you will see a message that says “connection successful.” You can then hit the button that says “lets go.”